

STEELCO DS500 DRS

Electrical connection available in 1x220V/60Hz or 3x208V/60Hz

[yourCEBA](#)

www.yourCEBA.com

STEELCO DS 500 DRS Single Door

Countertop - Undercounter - Surgical Instrument Washer Decontaminator



Product Support: 509-747-5027

[email:Support@washer-decontaminator.com](mailto:Support@washer-decontaminator.com)

Refer to websites for additional product information and literature.

[Website: yourCEBA.com](http://www.yourCEBA.com)

[Equipment WEBSITE: TBJ-SPD.com](#)

STEELCO DS 500 DRS Single Door

External dimensions (W x D x H) 900x630x840mm [35.43"x24.80"x33.07"] (with door closed)

Max rack dimensions (W x D x H) 520x500x530mm [20.47"x19.69"x20.84"]

DS500DRS in 3x208V,60Hz and UL listed		
STEELCO DS 500 DRS Single Door OPTIONS		
Flowmeter for chemicals		
Conductivity Sensor		
Glass window for DS500		
Stainless Steel Stand (900x560x600 mm) - 35.43" W x 22.04" D x 23.62" H - for Washer		
Stainless Steel Stand (600x560x600 mm) - 23.62" W x " x 22.04" D x 23.62" H - for Racks		
Additional Dosing Pump (3rd or 4th)		
Printer - Thermic, 40 columns, RS232 on front panel DS500		
DRAIN WATER only COOLING < 140° F (60° C)		
STEELCO DS 500 DRS Single Door RACKS - select at least 1 Rack		
1 level rack for lower level		C52
1 level rack for upper level with one wash arm underneath (can be added to the C52 Rack for 2 Level Rack)		C51
2 level rack for 4 instrument trays - Instrument Tray (full size) 19' x 9.85 x 2' (485x250x50 mm)		C98
3 level rack with wash arm for 6 instrument trays		C88
3 level with 10 Luer-Lock connections for cleaning cannulated instruments and working channels - 2 additional levels for trays		C48
4 level rack with 2 wash arms for 8 full size instrument trays		C78
MIS instrument rack 1 instrument tray 6 MIS instrument - 12 tubes & catheters w 10 Luer-Lock-connections		C95
Instrument Tray (full size) 19' x 9.85 x 2' (485x250x50 mm)		C62

STEELCO DS 500 DRS Standard Equipment

Pumps to Drain - does not require a dedicated floor drain - can drain into sink drain plumbing

Tank and air electrical heated

Manual operated hinged door

Re-circulating pump system

2 dosing pumps for chemical products

exhaust air condenser (steam condenser)

Cavity free wash chamber with rounded corners made of stainless steel AISI 316L

Drying unit of 150 mc/h [90cfm], equipped with double filter

Pre-filter F-class

HEPA Filter

Electrical air heater

Air drying for chamber and rack injectors

3-stage water filtering system

2 containers for chemical products.

40 Programs (20 pre-programmed)

STEELCO DS 500 DRS Single Door Dimensions: (W x D x H)

External dimensions 900x630x840mm [35.43"x24.80"x33.07"] (with door closed)

Max rack dimensions 520x500x530mm [20.47"x19.69"x20.84"]

STEELCO DS 500 DRS Electrical connection

available in 1x220V/60Hz or 3x208V/60Hz

Weight 250 pounds

Electrical tank heating 5 kW

Washing circuit with one pump for a total of 0.75 kW

Total power 5.6 kW

Water connection: cold, warm, DI

STEELCO DS 500DRS Single Door