

**STEELCO Washers**  
**Countertop Washers - Undercounter Washers**  
**Surgical Instrument Washers**  
**DS 500 DRS Specifications**

**STEELCO Countertop Washer Undercounter Washer DS500 DRS Surgical Instrument Washer Dimensions**

External Size Countertop Model: 35.43" W x 24.80" D x 33.07" H  
External Size Under Counter Model: 35.43" W x 24.80" D x 33.07" H  
Internal Rack Capacity: 20.47" W x 19.69" D x 20.87" H

**STEELCO Countertop Washer Undercounter DS500 DRS Surgical Instrument Washer Temperatures**

Cold Water Pre Wash Cycle: ambient to 95° F  
Wash Enzymatic/Detergent Cycle: 131° F to 137° F  
Elevated Thermal Rinse Cycle[s]: 194° F to 209° F  
Hot Air Dry Cycle: 158° F to 230° F

**STEELCO Countertop Washer Undercounter DS500 DRS Surgical Instrument Washer Cycle Times (plus Dry Time)**

Reservoir Capacity 3.5 gallons

20 pre-set programs

capacity for 20 additional customized programs

Rapid Wash: 8 minutes (with TOSI Test results "excellent")

Normal Wash: 12 to 15 minutes (range due to Temperature specified)

Intensive Wash: 22 minutes (reserved for heavy dried-on encrustation)

**STANDARD FEATURES for the Countertop Washer Undercounter DS500 DRS Surgical Instrument Washer:**

STEELCO Washer DS500 self cleaning chamber with rounded chamber corners

Temperature for washing and rinsing adjustable to 200 F

STEELCO Washer DS500 2 dosing pumps are standard additional are optional

LCD Display of processing time and remaining time

Complete Auto diagnosis for Trouble shooting system errors

Ss 232 connection for PC intemperance

STEELCO Washer DS500 integrated air exhaust steam condenser

Pump to Drain facilitates rapid cycle times

STEELCO Washer DS500 is UL listed

Available in 1x220V/60Hz or 3x208V/60Hz

STEELCO Washer DS500 HEPA ambient air Filter

**OPTIONS for the Countertop Washer Undercounter DS500 DRS Surgical Instrument Washer:**

STEELCO Washer DS500 Printer: cycle calculation and cycle tracking

Additional Dosing Pumps

STEELCO Washer DS500 Glass Window

Drain Water Cooling < 140° F (60° C)

1level rack for upper level with one wash arm underneath can be added to the C52 Rack for 2 Level Rack C51

1rack for lower level for STEELCO Washer DS500 C52

Instrument Tray (full size) 19' x 9.85 x 2' (485x250x50 mm) C62

2 level rack for 4 instrument trays for STEELCO Washer DS500 C98

3 level Micro Surgical Instrument and Cannulated Instrument Rack

with 10 Luer-Lock connections and 10 hose connections

on top level for cleaning lumens and working channels inside-and-out C48

3 level rack with wash arm for 6 instrument trays for STEELCO Washer DS500 C88

4 level rack with 2 wash arms for 8 full size instrument trays for STEELCO Washer DS500 C78

MIS instrument rack for one instrument tray with six MIS instrument for

with 12 tubes and catheter and 10 Luer-Lock-connections for cleaning lumens

and working channels inside-and-out C95

**STEELCO Washers**  
**Countertop Washers - Undercounter Washers**  
**Surgical Instrument Washers**  
**DS 500 DRS Specifications**