

STEELCO DS50, DS50DRS, DS50DRS-2

yourCEBA
www.yourCEBA.com

Countertop Washer Decontam **PASS-THRU or Single Door
 Surgical Instrument Washer - Washer Decontaminator - Washer Disinfector**

Product Support: 509-747-5027

[email:Support@washer-decontaminator.com](mailto:Support@washer-decontaminator.com)

Refer to websites for additional product information.

[Website: yourCEBA.com](http://www.yourCEBA.com)

[Equipment WEBSITE: TBJ-SPD.com](http://www.TBJ-SPD.com)



**STEELCO DS50, DS50DRS, and DS50DRS-2
 Chamber usable Size 15.75" x Wide 19.69" Deep 9.65 High**

STEELCO WASHER DS-50		
Countertop Tabletop Washer		
STEELCO WASHER DS-50-DRS		
Countertop Tabletop Washer with DRYER		
STEELCO WASHER DS-50-2 DRS 2 door Pass Through		
the only Countertop Tabletop Pass Through Washer with DRYER		
208 volt or 220 volt, 60Hz, UL Listed		
STEELCO DS50, DS50DRS, and DS50DRS-2 RACKS		
Vertical Instrument Basket	C03	
Small Tray Basket	C04	
Basic basket	C06	
Basket for small pieces	C28	
Insert for Instrument Containers with 3 positions for Cassettes 11 x 8 x 1.5	C90	
Insert for Instrument Containers with 4 positions for Cassettes 11 x 8 x 1.5	C91	
Tray for instruments	C92	

**STEELCO DS50, DS50DRS, and DS50DRS-2
 Pump to Drain - do not require a dedicated floor drain - pump into
 sink drain plumbing**

STEELCO DS50, DS50DRS, and DS50DRS-2 Dimensions: (W x D x H)

External dimensions 632x562x535mm
[24.9"x22.12"x21.06"] (with door closed)
 External dimensions 632x852x535mm
[24.9"x33.54"x21.06"] (with open door)
 Max rack dimensions 400x500x245mm
[15.75"x19.69"x9.65"]

Max. basket capacity 49 lt. [13gal.]
Electrical connection 1x208V/60Hz
Electrical tank heating 3 kW

Washing circuit with one pump for a total of 110W
Total power 3.5 kW
Water connections: one

STEELCO DS50, DS50DRS, and DS50DRS-2 Description:

washing chamber entirely made of stainless steel AISI 316L
smooth chamber design to eliminate dead legs for potential bacteria grow
self cleaning chamber with rounded corners
water filtering system for excellent washing results
Pass through version with two manual bottom hinged doors
Frontal basket loading
Automatic Rack coupling system for injector cleaning and drying to guarantee
highest quality cleaning and efficient drying
cleaning and drying results

- Full drainage of re-circulating pump
- Thermal decontamination possibilities up to 93°C [200°F]
- PT1000 temperature probe
- Equipped with 2 dosing pumps for chemicals, microprocessor controlled.
- Level control for chemical jugs with minimum level alarm
- LED display as user interface
- 20 Pre-set programs

Total control of cycle parameters:
time, temperature, chemical dosing, water quantity,
drying

- Exhaust air condenser
- Equipped with hot air dryer

not included in price:

- > start-up and training
- > installation
- > transport & insurance