

HYDRIM Cleaning Solution

Enzymatic Surgical Instrument Cleaning Detergent with Lubricant

The UPS ergo-2 is for use in HYDRIM Washers, such as the L110w and M2

The all-in-ONE HYDRIM Cleaning Solution Enzymatic Detergent Cleaner and Lubricant delivers:


- multi-tiered enzymatic detergent complex
- stain removers for residue free cleaning,
- surgical instrument lubricant
- 100% biodegradable and phosphate free,
- neutral pH enhances surgical instrument stainless steel passive layer
- water softeners remove mineral encrustations

The HYDRIM Cleaning Solution Enzymatic Surgical Instrument Detergents **cut costs**.

The Surgical Instrument Detergent Enzyme Lubricants clean **residue free**.

Four Enzyme Surgical Instrument Cleaning Detergent Lubricants **clean faster**.

- Enzyme Cleaners (protease, amylase, carbohydrase, and lipase)
- Detergent Cleaners (surfactant cleaners)
- Lubricants (water soluble and free rinsing)
- Surface Conditioners (strengthen the chromium passive layer of stainless steel)



The all-in-ONE HYDRIM Cleaning Solution Enzyme Detergent Lubricant Cleaners delivers the:

- highest concentration enzymatic detergent available,
- combined enzyme and detergent cleaning treatments that improve cleaning outcomes,
- increased performance of surgical instruments and an extended functional lifespan.

The Neutral pH (7) Surgical Instrument Enzymatic Detergent Lubricants clean fast and residue free.

The all-in-ONE Surgical Instrument Cleaners deliver the four enzymes that are necessary to break down all forms of proteinaceous bioburden and the surfactant detergents to clean the surface. The all-in-ONE is highly concentrated for faster cleaning and lower cleaning costs.



Customer Service and Purchasing: 1-509-981-1553 email: Info@yourCEBA.com
Product Technical Support: John Temple, 1-509-747-5027, email: Temple@yourCEBA.com
Websites: www.yourCEBA.com