

HYDRIM Cleaning Solution

Enzymatic Surgical Instrument Cleaning Detergent Lubricant

UPS ergo-2 for use in HYDRIM Washers (i.e.: L110w and M2)


The all-in-ONE HYDRIM Cleaning Solution Enzymatic Detergent Cleaner Lubricant delivers:

- multi-tiered enzymatic detergent complex
- Enzyme Cleaners (protease, amylase, carbohydrase, and lipase)
- Detergent Cleaners (surfactant cleaners)
- Lubricants (water soluble and free rinsing)
- Surface Conditioners (strengthen the chromium passive layer of stainless steel)
- stain removers for residue free cleaning,
- neutral pH enhances surgical instrument stainless steel passive layer
- water softeners remove mineral encrustations

The HYDRIM Cleaning Solution Enzymatic Surgical Instrument Detergents **cut costs**.

The Surgical Instrument Detergent Enzyme Lubricants clean **residue free**.

Four Enzyme Surgical Instrument Cleaning Detergent Lubricants **clean faster**.



The Neutral pH (7) HYDRIM Cleaning Solution clean fast and residue free. The all-in-ONE HYDRIM Cleaning Solution delivers the four enzymes necessary to break down all forms of proteinaceous bioburden and the surfactant detergents to clean the surface. The all-in-ONE HYDRIM Cleaning Solution is highly concentrated for faster cleaning and lower cleaning costs.



Customer Service: 1-509-220-2456 CustomerService@ergo-Logistics.com
Admin and Technical Support: 1-509-747-5027 Temple@ergo-Logistics.com
Website: www.ergo-Logistics.com